

**Mercaptan sulphur in
Gasoline and Kerosene
Aviation Turbine and
Distillate Fuels according to
ASTM D 3227**

Application

Use

Determination of mercaptan sulphur in a range from 3 mg/kg – 100 mg/kg.

Appliances

- Titrator: TL 7000/TL 7750 M1
- Basic device
- Magnetic stirrer TM 235
- 10 mL Exchange unit WA 10, with amber glass bottle for the titrant, complete

Electrodes

- Electrode: AgS 62 RG or AG 1100 (sulphidized) + A 1180/H 1180 glass electrode
- Electrode cable: L 1 A or L 1 A + L 1 N

Reagents

- Titrant: AgNO₃ 0.01 mol/l in isopropanol (IPA)
- Titer determination: with KI or NaCl solution 0.1 mol/l
- Solvent: Sodium acetate trihydrate solution in IPA

Description

Preparation and standardization of the K-I and alcoholic AgNO₃ solutions

Dissolve 17 g (weigh to 0.01g) of KI in 100 ml of water in a 1 L volumetric flask and dilute to 1 L. Calculate the exact molarity. It is also possible to use NaCl instead of KI. Commercial available 0.1 mol/L solutions for NaCl can be also used.

Use a standard ampoule **0.1 m AgNO₃** and fill it up with IPA (99%) in a 1 L volumetric flask.

Add 6 drops of conc. HNO₃ to 100 ml of water in a 250 ml tall form beaker. Remove oxides of nitrogen by boiling for 5 min. Cool to ambient temperature. Pipette 5 ml of a 0.1 m KI solution into the beaker and titrate with the 0.1 m AgNO₃ to an inflection point.


The **0.01 m AgNO₃** should be prepared daily by dilution of the 0.1 m standard. Calculate the exact molarity.

Application

Method for titer determination

Page 1

GLP documentation



Method data

Method name:	AgNO3 0,01 M in IPA	Titration duration:	7 m 7 s
End date:	25.03.15	End time:	17:31:48

Titration data

Sample ID:	NaCl 0,4082 g/100 ml	Pattern:	1.000 ml
Start mV:	-46,6 mV	End mV:	-267,9 mV
EQ:	6.970 ml / -221,1 mV	Titer:	0.01002 mol/l
Mean value:	---	RSD:	---

Calculation formula

Titer:	$(V^2 F_2) / ((EQ_1 - B)^2 M^2 F_1) \rightarrow WA$	Mol (M):	1.00000
Pattern (V):	1.000 ml (m)	Factor 2 (F2):	69.8458
Blank value (B):	0.0000 ml	Factor 1 (F1):	1000.0000
Statistics:	1 from 3		

Application

Page 2

Method data overall view

Method name:	AgNO3 0,01 M in IPA	Created at:	03/25/15 17:24:34
Method type:	Automatic titration	Last modification:	03/25/15 17:24:34
Measured value:	mV	Damping settings:	strong
Titration mode:	Dynamic	Documentation:	GLP
Dynamic:	User-defined:	Max. step size:	0.5000 ml
		Slope max ml:	15.00 mV/min
		Min. step size:	0.0200 ml
		Slope min ml:	230.00 mV/min
Measuring speed / drift:	User-defined:	minimum holding time:	05 s
		maximum holding time:	15 s
		Measuring time:	04 s
		Drift:	05 mV/min
Initial waiting time:	0 s		
Titration direction:	Decrease		
Pretitration:	Off		
End value:	Off		
EQ:	On (1)		
Slope value:	User-defined	Value:	400

Dosing parameter

Dosing speed:	100.00 %	Filling speed:	30 s
Maximum dosing volume:	10.00 ml		

Unit values

Unit size:	10ml
Unit ID:	10035433
Reagent:	AgNO3 in IPA
Batch ID:	no entry
Concentration [mol/l]:	0.01000
Determined at:	03/25/15 22:33:43
Expire date:	--
Opened/compounded:	--
Test according ISO 8655:	05/04/12
Last modification:	03/25/15 15:33:54

Device information

Application

Preparation of the Solvent

Dissolve 2.7 g of sodium acetate trihydrate in 20 ml oxygen free water and pour into 975 ml 2-propanol (IPA). Add 4.6 ml of glacial acetic acid. Remove dissolved oxygen with a rapid stream of nitrogen for 10 min each day prior to use. Keep protected from the atmosphere.

Preparation of the CdSO₄ solution

Dissolve 150 g of CdSO₄ (3CdSO₄ * 8 H₂O) in water. Add 10 ml of H₂SO₄ (1:5) and dilute to 1L with water

Connection of the electrode

The AgS 62 RG is directly connected to pH/mV socket with cable L 1 A.

Option:

The A 1180 is connected with the cable L 1 A to pH/mV socket". The Ag 1100 (sulphidized) is connected with cable L 1 N to the reference socket.

Titration

Preparation of the sample: Removal of Hydrogen Sulfide

Test the sample quality for Hydrogen sulfide (H₂S) by shaking 5 ml of the sample with 5 ml of the acid CdSO₄ solution. If no precipitate appears, proceed with the analysis of the sample described below. If a yellow precipitate appears, remove the H₂S in the following matter: Place a quantity of the sample (i.e 200 ml) , three, two, four times that required for the analysis, in a separatory funnel containing a volume of the acid CdSO₄ solution equal to one half of the sample (i.e. 100 ml) and shake vigorously. Draw off and discard the aqueous phase, and wash the sample with three 25-30 ml portions of water, withdrawing the water after each washing. Repeat the extraction with CdSO₄ until all of the H₂S has been removed.


Measure with a pipet or weigh 20 to 30 ml of the original or treated sample into a 150 ml titration beaker (titration vessel) containing 70 ml of the solvent mixture. Immediately immerse the electrodes and buret tip into the sample. Titrate with the 0.01 m AgNO₃ with the attached titration parameters. After the titration the electrodes should be rinsed with alcohol and then with water.

Application

Result page 1

GLP documentation

Titration graph



Method data

Method name:	R-SH without H2S	Titration duration:	6 m 40 s
End date:	25.03.15	End time:	19:09:14

Titration data

Sample ID:	Petroleum Standard	Pattern:	20.000 ml
Start mV:	400.7 mV	End mV:	-213.9 mV

EQ:	1.121 ml / -32.0 mV	R-SH:	21.9 ppm
-----	---------------------	-------	----------

Calculation formula

R-SH:	$(EQ1-B)^T * M * F1 / (V * F2)$	Mol (M):	32.06000
Blank value (B):	0.0000 ml	Titre (T):	0.00998 (a)
Factor 1 (F1):	1000.0000	Pattern (V):	20.000 ml (m)
Factor 2 (F2):	0.8200	Statistics:	Off

Application

Result page 2:

Method data overall view

Method name:	R-SH without H2S	Created at:	03/25/15 19:02:24
Method type:	Automatic titration	Last modification:	03/25/15 19:02:24
Measured value:	mV	Damping settings:	strong
Titration mode:	Linear	Documentation:	GLP
Linear steps:	0.050 ml		

Measuring speed / drift:	User-defined:	minimum holding time:	05 s
		maximum holding time:	15 s
		Measuring time:	04 s
		Drift:	05 mV/min
Initial waiting time:	0 s		
Titration direction:	Decrease		
Pretitration:	Off		
End value:	Off		
EQ:	On (1)		
Slope value:	User-defined	Value:	500

Dosing parameter

Dosing speed:	100.00 % (20.00 ml/min)	Filling speed:	30 s
Maximum dosing volume:	10.00 ml		

Unit values

Unit size:	10ml
Unit ID:	10035433
Reagent:	AgNO3 in IPA
Batch ID:	no entry
Concentration [mol/l]:	0.00998
Determined at:	03/25/15 18:22:34
Expire date:	--
Opened/compounded:	--
Test according ISO 8655:	05/04/12
Last modification:	03/25/15 18:22:40

Application

Notes

If you have any questions on the application, you can feel free to contact us..

SI Analytics GmbH
Hattenbergstr. 10
55122 Mainz
Germany

Phone: +49 (0) 6131 / 66 – 5062
+49 (0) 6131 / 66 – 5118
Fax: +49 (0) 6131 / 66 – 5001
E-Mail: titration@si-analytics.com
Homepage: www.si-analytics.com