

## Flow Solution™ FS 3700 Automated Chemistry Analyzer


Methylene Blue Active Substances (MBAS) Index using Continuous Flow Analysis  
ISO 16265  
Cartridge Part Number 330358CT

### Scope and Application

This method is used for the determination of the methylene blue active substances (MBAS) index in drinking water, ground water, surface water, domestic and industrial wastes according to **ISO Method 16265**. Anionic surfactants are the most important substances showing methylene blue activity. This method is, therefore, useful for estimating the anionic surfactant content (e.g. soaps) of water.

### Method Performance

<b>Range</b>	0.025 - 2.0 mg/L as LAS
<b>Rate</b>	24 samples/hour
<b>Precision</b>	<5% RSD at mid-point of range
<b>Method Detection Limit (MDL)</b>	0.008 mg/L as LAS



**Figure 1.** General Flow Diagram for Methylene Blue Active Substances by ISO 16265

## Reagents and Calibrants

Chemical Name	CAS #	Chemical Formula	Part Number
Boric acid	10043-35-3	$H_3BO_4$	
Brij®-35	9002-92-0	$(C_{12}H_{26}O)_n$	326126
Chloroform	67-66-3	$CHCl_3$	
Ethanol	64-17-5	$C_2H_6O$	
Hydrochloric acid	7647-01-0	HCl	
Isopropyl alcohol	67-63-0	$C_3H_8O$	
Methylene blue	7220-79-3	$C_{16}H_{18}N_3S$	
Petroleum ether	8032-32-4		
Sodium dodecylbenzenesulfonate	25155-30-0	$CH_3(CH_2)_{11}C_6H_4SO_3Na$	
Sodium dodecyl sulfate	151-21-3	$CH_3(CH_2)_{11}OSO_3Na$	
Sodium hydroxide	1310-73-2	NaOH	
Sodium phosphate, monobasic	7558-80-7	$NaH_2PO_4$	
Sulfuric acid, concentrated	7664-93-9	$H_2SO_4$	
Water, deionized		$H_2O$	
<b>Additionally, the following chemical may be needed for sample preservation or treatment.</b>			
Formaldehyde	50-00-0	$CH_2O$	

## Summary of ISO 16265

### Method

- When such anionic surfactants are mixed with an alkaline methylene blue solution, an ion pair is formed. In mixtures of chloroform and water, the ion pairs will extract into the chloroform layer, transferring the blue color into the organic phase (the unassociated dye has an extremely small solubility in chloroform). Separating the organic phase from the aqueous phase removes interferences and yields chloroform containing the extracted methylene blue-anionic complexes. The organic phase is then treated with an acidic methylene blue solution and the absorbance measured at 640 nm, which is the absorbance maximum of methylene blue in chloroform.
- The anionic surfactants most often found in wastewater are soluble sodium salts of the alkyl sulfates and the alkyl benzene sulfonates.
- Separation of the aqueous and organic phases can be accomplished using different methods including the use of hydrophobic membranes. For proper operation of this method, ensure that only the organic phase reaches the detector.

## Interferences

1. Anionic surfactants are common components of detergents and may be present in residual amounts on laboratory glassware. Wash all glassware with dilute acid prior to usage to prevent interference from contamination.
2. Methylene blue solutions must be protected from oxidizing substances.
3. Filtration of samples containing particulates may be required. As a result, some anionic surfactant loss may occur.
4. Heavy metals, cationic surfactants, or other compounds (e.g. amines) that can compete with methylene blue to form an ion pair with MBAS will reduce response.
5. High concentrations of chloride, bromide, and nitrate interfere by forming ion pairs with methylene blue. These interferences are removed using a second extraction at an acidic pH.<sup>2</sup>
6. Other compatible procedures for removing or suppressing interferences may be used, provided they do not adversely affect overall method performance.




Figure 2. MBAS Calibration Series




Figure 3. Calibration Curve and Statistics

Table 2. Consumables and Spare Parts for Methylene Blue Active Substances, ISO 16265

Consumables	Part Number
Pump tubes kit - Methylene Blue Active Substances, ISO 16265	330358TK
Viton® pump tube assembly – White/White	319714
Viton® pump tube assembly – Grey/Grey	319713
Sample Vials, Glass 8-mL (13 X 100 mm) (pack of 1000)	A000514
Phase Separation Membrane (5 pk)	A001520
PEEK Autosampler Probe for RA/3090/3360 Sampler	325331
Brij®-35	326126

Optional Accessories	Part Number
Teflon® Phase Separation Membrane (25 pk)	A002040
Plastic Membrane Tweezers – Wide, Flat Tip	328929
Debubbler, Polysulfone for the FS3100, FS IV, FS 3700	A000172
FS3700 8-port FIA Valve & Cable	330394
400 µL Sample Loop/Bypass Loop for 8-port valve	319334
Pump Tube Assortment Pack (2 of every-sized pump tube)	A000362

Normal pump tubes should be replaced monthly, Viton pump tubes may need to be replaced weekly or as-needed to maintain system performance. Maximum life expectancy for normal pump tubes is approximately 800 hours.


  
 151 Graham Road  
 PO Box 9010  
 College Station, Texas  
 77842-9010  
 (979) 690-1711  
 (800) 653-1711 USA/Canada  
 (979) 690-0440 Fax  
 www.oico.com  
 E-mail: OI-Mail@Xylem.com